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
930 Compact IC Flex

Compact ion chromatography system for
routine analysis

 **Metrohm**

Handwritten signature

Handwritten signature

930 Compact IC Flex: Entry-level model and workhorse for routine analysis


02

The 930 Compact IC Flex is Metrohm's new ion chromatograph for routine analysis developed with a focus on the requirements of contract laboratories and QC laboratories in all kinds of industries. Rugged design, perfect ease of use and outstanding system reliability are key features of the 930 Compact IC Flex. Monitoring and control functions for system parameters, service intervals, calibrations, and results take the stress out of daily routine operation while ensuring high-quality measuring results – even when things get hectic in the laboratory.


The new 930 Compact IC Flex system – as its name says – offers the highest possible degree of flexibility: You can choose from a complete range of separation columns, suppressors and detectors to configure a customized so-

lution that meets your specific analytical requirements. If necessary, the 930 Compact IC Flex can also be fitted with a column oven as well as an eluent and sample degasser. It goes without saying that the unique Metrohm Inline Sample Preparation («MISP») techniques are available in the new system, as are numerous options for automation.

Your way to your customized 930 Compact IC Flex system is simple: Use our new online configurator (ic930.metrohm.com) to select from a wide range of options and put together precisely the right system to meet your requirements. With the new 930 Compact IC Flex the dream of a customized high precision tool for routine analysis has at last come true!

The new 930 Compact IC Flex configurator allows you to put together your customized IC system for routine analysis in just a few steps. Combine up to 90 different instruments and accessory parts depending on the requirements of your application. Try it out on ic930.metrohm.com

03

Highlights

- Compact system for routine analysis of anions, cations and polar substances in the range of $\mu\text{g/L}$ to g/L
- Intelligent ion chromatography for superior reliability
- Modular kit of system components and accessories for custom system configuration
- Can be combined with all types of detection: conductivity, UV/VIS, amperometry
- Space-saving design, easily accessible system components
- STREAM – the green way of suppression
- Complete automation and unique Metrohm Inline Sample Preparation («MISP») possible
- Complies with all GLP and FDA requirements
- Multi-language MagIC Net software for simple and intuitive operation
- Comprehensive monitoring- and control functions for high quality results

Handwritten signature

Handwritten signature

Applications


04

Small footprint, competitive price and measuring results of outstanding quality – Metrohm has revolutionized ion chromatography with the introduction of compact IC instruments. Nowadays these instruments dominate the field in routine water and environmental analysis. They are used for the investigation of drinking, surface, ground and waste waters.

However, the new 930 Compact IC Flex was not developed solely for these industries. Thanks to its great reliability, the system is perfectly suited for use in the pharmaceutical industry, e.g. for analyzing infusion solutions. Thanks to the system's flexibility, the 930 Compact IC Flex instruments are however also suitable for the chemical and food industries, where (in addition to conductiv-

ity detection) amperometric and UV/VIS detection are regularly used. Small to medium-sized laboratories on the other hand benefit particularly from the 930 Compact IC Flex' excellent price-performance ratio and the system's ease of use.

Furthermore, low detection limits also make the instruments in the 930 Compact IC Flex family an excellent choice for routine analysis in power plants with detection limits down to trace levels. And finally, the new compact ion chromatographs from Metrohm are not only suitable for routine analysis in the petrochemical industry but also for the quality monitoring of alternative fuels, e.g. bio-ethanol and biodiesel.

The 930 Compact IC Flex can be used to analyze gaseous, liquid, and solid samples. The Combustion IC system shown can be used for differentiated determination of the halogens and sulfur in combustible samples, e.g. plastics, raw or end products in the petroleum industry, samples from waste management or electronic components.

Fully automated analyses for more work efficiency and analysis reliability

Automation saves time and money

The 930 Compact IC Flex offers completely automated operation. For liquid samples alone, there are six different autosamplers available in a total of 23 versions. The individual systems differ with respect to sample capacity, cooling, liquid handling functions and additional valve options. This means it is always possible to find the optimum automation tailored to meet your requirements.

Intelligent injection techniques for an extended working range





Apart from the full-loop and internal-loop injection, the 930 Compact IC Flex can also be combined with various intelligent injection techniques such as the «MiPT» (Metrohm intelligent Partial-Loop technique) and «MiPuT» (Metrohm intelligent Pick-up technique). The variable injection volume of MiPT covers a sample measuring range that extends across 4 orders of magnitude. This means that samples in the range of 10 µg/L to 100 mg/L can be analyzed with a single calibration.

Metrohm Inline Sample Preparation («MISP»): More efficiency in routine operation and new fields of application

The Metrohm Inline Sample Preparation techniques – some of which are patented – significantly expand the scope of application for ion chromatography. Inline Ultrafiltration, Inline Dilution, Inline Dialysis or other techniques make even the most challenging samples manageable, such as suspensions or waste water samples that are loaded with proteins or extremely contaminated. The combination of Inline Dilution and Inline Ultrafiltration stands out here in particular as one of the most frequently used routine applications.

Anion and cation determination with just one autosampler

The 930 Compact IC Flex system allows simultaneous analysis of anions and cations down to the µg/L range. A setup of this kind is comprised of two Compact IC instruments sharing an autosampler. Together, they form a fully automated analysis system for processing sample series 24/7 determining a wide range of ionic components.

930 Compact IC Flex – the ideal partner for routine analysis

06

Maximum reliability

The 930 Compact IC Flex excels by superior reliability. The system is self-monitoring, i.e.

- System components are immediately recognized
- Instruments and other parameters are automatically integrated into the method
- All system and method parameters are monitored permanently
- Measuring results are traceable to every single step of the analysis

If a parameter exceeds a defined limit, the system automatically tells you so sending a message – in plain text. Operator errors that could theoretically cause damage to the separation column, for example, are thus virtually ruled out in practice


STREAM (Suppressor Treatment with Reused Eluent After Measuring) – the green way of suppression

Whether sequential, chemical or without any suppression: You have freedom of choice with the 930 Compact IC Flex. The suitable rotor is selected depending on the application and the column dimensions: «MSM-HC», «MSM II» or «MSM-LC». Each of these rotors fits in the transparent suppressor housing and, due to the robust design, comes with a 10-year manufacturer's warranty.

All 930 Compact IC Flex versions with suppression are equipped with STREAM. In the STREAM setup, the suppressed eluent is used for rinsing the regenerated suppressor unit after the detection. This means no additional rinsing medium is necessary. Apart from that, the flow of regenerant can be reduced to a minimum. The benefits are less need of chemicals and less liquid waste. Furthermore, STREAM enables nonstop system operation for at least two weeks – or even longer than that with inline preparation of the regenerant. This saves on manual working steps, reduces maintenance and thus helps cutting running costs.


Professional Liquid Handling

A peristaltic pump and the patented 800 Dosino are available for transporting auxiliary solutions in sample preparation, for transferring samples and for rinsing or regeneration procedures. The 930 Compact IC Flex can manage entirely without a peristaltic pump for suppressor regeneration. In this case, we recommend the new Dosino Regeneration «DR». This reduces the need for maintenance and increases system reliability.

Working continuously without manual intervention

Combined with the 941 Eluent Production Module, the 930 Compact IC Flex integrates automated inline preparation of eluents of any composition and concentration. Connect an ultrapure water system (e.g.; ELGA PURELAB flexS/6) to the system and you may use conventional tap water for your ion chromatography. Automated inline eluent preparation ensures stable retention times, contamination-free working and saves manual working steps.

07

Perfect flexibility for the best application solution

Each sample may require different sample preparation and/or analysis. This is why the 930 Compact IC Flex offers a wide range of system components that can be selected to meet the particular purpose and requirements: with or without column oven; Dose-in Gradient; eluent and sample degasser; conductivity, UV/VIS or amperometric detection or a kind of inline sample preparation it may be – the optimum solution is available with the 930 Compact IC Flex. Moreover the 930 Compact IC Flex can be operated with any separation column, regardless of the base material, particle size or dimensions.

MagIC Net – user-friendly software

The instruments of the 930 Compact IC Flex series are controlled by the proven MagIC Net ion chromatography software. Freely configurable user windows, and graphic symbols for the individual system components make the software simple and intuitive to use. MagIC Net is available in no fewer than 14 languages!

Reliable results – automatically!

A wide range of monitoring and control functions ensure highest reliability – of both the system itself and the quality of the results produced. Be it the number of injections on a separation column, a parameter exceeding defined tolerance limits for results or checking the calibration with a check standard – the 930 Compact IC Flex provides complete information. That's not all: If required, the system intervenes and takes action automatically. For example, recalibration is carried out automatically if the check standard should fall outside the defined limits.

MagIC Net – the ion chromatography software


08

The instruments in the 930 Compact IC Flex series are controlled by the proven MagIC Net ion chromatography software. MagIC Net also controls any peripheral devices for Liquid Handling and automation. MagIC Net records the results produced, enables modern data management and reporting as required by the user.


MagIC Net is easy to command. The user interface can be freely configured and adapted to the needs of the user. Thus, only those windows are visible that are actually needed; the kind and scope of information in these windows can, in turn, be defined as required by the user. If required, system command can be simplified to a single click on the start or stop button on the screen! As MagIC Net is available in 14 languages, linguistic misunderstandings and resulting errors by the operator are virtually ruled out.

MagIC Net provides self-monitoring of the system and ensures that any results produced are checked automatically. Thus, the software makes logical decisions and takes action by itself, if required. A good example would be the determination of the optimum dilution factor: If the concentration of the analyte is outside the calibrated range, then the system automatically calculates the required dilution factor and initiates dilution of the sample, ensuring that any results produced are always reliable.

Mag Net

Technical information

09

General	• Compact IC system with modular design • Intelligent system components • Combination with various detectors possible • Metal-free flow path • 3-year warranty
Intelligent system components	Intelligent technology is integrated in the following system components, among others:
Eluent and sample degasser	• iPump • iDetector • iColumn • 800 Dosino
High-pressure pump	Organic modifier 0 - 100% (no PFC [perfluorocarbons]) Material fluoropolymer Serial dual-piston pump with two valves and flow range-optimized, intelligent pump heads Flow rate 0.001 - 20 mL/min
Injection valve	Injection volume: Internal loop 0.25, 1 µL Sample loops 1.5, 5, 10, 20, 50, 100, 250, 1'000 µL
Column oven	Temperature range 0...+80 °C (ambient temperature +5...+40 °C) Stability <0.05 °C
Suppressors	«MSM II», «MSM-HC» and «MSM-LC» Metrohm Suppressor Modules for chemical suppression Type «Micro Packed Bed» suppressor Regeneration STREAM with peristaltic pump or 800 Dosino Organic modifier 0 - 100% Warranty 10 years «MCS» Metrohm CO₂ Suppressor Type CO₂ removal with fluoropolymer technology Organic modifier 0 - 100% (no PFC [perfluorocarbons])
Detectors	Options for integration in the system include: • Conductivity detection • UV/VIS detection • Amperometric detection
Conductivity detector	Intelligent high-performance conductivity detector with DSP – «Digital Signal Processing» Measuring range 0 - 15'000 µS/cm – without range switching Temperature 20 - 50 °C in 5 °C increments Temperature constancy <0.001 °C Cell volume 0.8 µL Electronic noise <0.1 nS/cm (at 1 µS/cm) Baseline noise <0.2 nS/cm (e.g. A Supp 5, standard conditions)
Peristaltic pumps	Rotational speed 0 - 42 rpm in increments of 6 rpm Shift direction clockwise and counterclockwise
Gradients	Dose-in Gradient: binary, ternary, quaternary, quinary Progression: step, linear
Automation	Combinable Autosamplers: 858 Professional Sample Processor, 919 IC Autosampler plus, 863 Compact IC Autosampler, 889 IC Sample Center, 814 USB Sample Processor, 815 Robotic USB Sample Processor XL
Metrohm Inline Sample Preparation «MISP»	Options for integration in the system include: • Inline Ultrafiltration • Inline Dialysis • Inline Matrix Elimination • Inline Dilution • Inline Extraction
Metrohm injection techniques	Every 930 Compact IC Flex is pre-installed with full-loop injection; options for integration in the system include: • Internal-loop injection • Metrohm intelligent Partial Loop Injection «MIPT» • Metrohm intelligent Pick-up injection «MiPuT»
Control	MagIC Net Compact, Professional and Multi ion chromatography software

Handwritten signature

Handwritten signature

Ordering Information

10

	930 Compact IC Flex instruments
2.930.1100	930 Compact IC Flex
2.930.1160	930 Compact IC Flex Deg
2.930.1200	930 Compact IC Flex ChS
2.930.1260	930 Compact IC Flex ChS/Deg
2.930.1300	930 Compact IC Flex ChS/PP
2.930.1360	930 Compact IC Flex ChS/PP/Deg
2.930.1400	930 Compact IC Flex SeS
2.930.1460	930 Compact IC Flex SeS/Deg
2.930.1500	930 Compact IC Flex SeS/PP
2.930.1560	930 Compact IC Flex SeS/PP/Deg
2.930.2100	930 Compact IC Flex Oven
2.930.2160	930 Compact IC Flex Oven/Deg
2.930.2200	930 Compact IC Flex Oven ChS
2.930.2260	930 Compact IC Flex Oven/ChS/Deg
2.930.2300	930 Compact IC Flex Oven ChS/PP
2.930.2360	930 Compact IC Flex Oven ChS/PP/Deg
2.930.2400	930 Compact IC Flex Oven SeS
2.930.2460	930 Compact IC Flex Oven/SeS/Deg
2.930.2500	930 Compact IC Flex Oven SeS/PP
2.930.2560	930 Compact IC Flex Oven SeS/PP/Deg

Detection



2.850.9010	IC Conductivity Detector
2.850.9110	IC Amperometric Detector
2.944.0010	944 Professional UV/VIS Detector Vario
2.945.0010	945 Professional Detector Vario – Conductivity
2.945.0020	945 Professional Detector Vario – Amperometry
2.945.0030	945 Professional Detector Vario – Conductivity & Amperometry

MagIC Net software

6.6059.301	MagIC Net 3.0 Compact
6.6059.302	MagIC Net 3.0 Professional
6.6059.303	MagIC Net 3.0 Multi

Automation

2.814.0130	814 USB Sample Processor – 2T
2.815.0130	815 Robotic USB Sample Processor XL – 2T
2.858.0010	858 Professional Sample Processor
2.858.0020	858 Professional Sample Processor – Pump
2.858.0030	858 Professional Sample Processor – Pump – Injector
2.863.0010	863 Compact IC Autosampler
2.889.0010	889 IC Sample Center
2.889.0020	889 IC Sample Center – cool
2.919.0020	919 IC Autosampler plus


858 Professional Sample Processor – selection of accessories

6.2041.440	Standard rack 148 × 11 mL + 3 × 300 mL
6.2041.760	Rack 54 × 11 mL + 1 × 300 mL
6.2041.480	Rack 159 × 2 mL + 3 × 300 mL
6.2743.050	PP sample vessels (11 mL); 2'000 units
6.2743.040	PP sample vessels (2.5 mL); 2'000 units
6.2743.070	PP stopper with perforation, for sealing the sample vessels; 2'000 units
6.5330.130	IC Equipment: Liquid Handling Station

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AMM

Automation in ion chromatography

Save time and money through automated sample
preparation and analysis

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Ω Metrohm




Blas



Automation for more analytical reliability

02

In ion chromatography anions, cations, carbohydrates and polar substances are determined in a variety of matrices. Depending on the application, solid, liquid or gaseous samples need to be analyzed. Complete automation of such determinations is increasing in importance all the time. For good reasons: Automation reduces manual steps to a minimum providing for better reproducibility and accuracy of results.

Metrohm's modern automation solutions for ion chromatography are highly flexible. Sample volumes ranging from one microliter to half a liter can be handled without any problem. Even special sample vessels can be placed directly on the sample changer, as we make our sample racks according to our customers' individual requirements. This eliminates the risk of contamination that occurs when the sample is transferred into another vessel. The combination of sample introduction with the different Metrohm Inline Sample Preparation and intelligent injection techniques saves time and cost and makes the entire analysis, including sample preparation, traceable right down to every single step.

01/11/14

Highlights

- Fully automatic sample introduction
- High precision and accuracy of results
- Cost and time savings
- Parallel working
- Combination with Metrohm Inline Sample Preparation techniques
- Professionell Liquid Handling and intelligent injection techniques
- Flexibility in the volume of samples
- Use of any type and size of sample vessels
- Completely metal-free sample channel
- Optional cooling function
- Sample rack and possibility for rinsing the sample needle can be selected freely
- Each step is freely programmable with the ion chromatography software MagIC Net
- Robust, reliable and continuous handling of complete sample series

01/11/14

Metrohm Inline Sample Preparation techniques

04

Crystal clear, clinically pure, free of particles – this is the description of the ideal sample for ion chromatography. But things are usually very different in real life. Various sample preparation steps are necessary in order to protect the separation column and the entire IC system from precipitation and blockages.


Metrohm enables carrying out these sample preparation steps inline in the course of the flow path, i.e. already during sample introduction. Unique techniques – some of them patented by us – are available for this. These techniques reduce manual steps, are very robust, reliable and require hardly any maintenance.

Metrohm Inline Sample Preparation opens up new fields of application. Thus, the reliable elimination of potential

sources of contamination opens up the field of ultratrace analysis. Moreover, each single step is traceable. Users save valuable time and increase analytical reliability.

Methods of Metrohm Inline Sample Preparation (MISP)

- Inline Ultrafiltration
- Inline Dialysis
- Inline Dilution
- Inline Extraction
- Inline Matrix Elimination
- Inline Neutralization
- Inline Cation Removal
- Inline Preconcentration
- Inline Calibration
- Inline Spiking

The dialysis cell shown is the heart of Metrohm Inline Dialysis. This patented technique separates not only particles from your analytes, but also colloids, oil components and large molecules, such as proteins. Thus you can determine anions, cations, carbohydrates and small polar substances in the $\mu\text{g/L}$ to g/L range in the following sample matrices:

- emulsions and dispersions
- process water, washing water and waste water
- oil-containing samples such as cutting oils and petroleum-containing samples
- samples from fermentation processes
- dairy products and other protein-containing samples
- body fluids
- organically contaminated samples
- concentrated fruit and vegetable juices


Professional Liquid Handling guarantees correct results

A scrupulous, contamination-free working method is fundamental for any Liquid Handling in ion chromatography. Moreover, particular attention must be paid to accuracy and precision. In order to achieve this, Metrohm invented the patented 800 Dosino. This dosing system can aspirate liquids fast and dose them with extreme precision (to an accuracy of 0.2 µL). Thus, the 800 Dosino guarantees flow and dosing accuracy over a long period. Besides the classical full-loop techniques, alternative injection techniques are also feasible:


- Internal Loop Injection Technique
- intelligent Partial Loop Injection Technique (MIPT)
- intelligent Pick-up Injection Technique (MIPuT)

For ultratrace analysis the Dosino is combined with a buffer volume and a trap column. This sophisticated method prevents carryover from sample to sample and also prevents effects of auxiliary solutions on the chromatogram. As a result, ions can be determined reliably down to the single-digit ng/L range.

05


IC system for determining ultratraces in strongly alkaline or acidic samples: The IC system is equipped with Inline Neutralization and Inline Preconcentration. Inline Neutralization ensures that even strongly alkaline samples can be preconcentrated and thus replaces time-consuming and cost-intensive manual sample preparation with solid-phase extraction cartridges (SPE cartridges). Variable Inline Preconcentration enables reliable calibration with a single standard. By combining both techniques with the professional Liquid Handling of the 800 Dosino it is possible – amongst other things – to reliably analyze ultratraces of anions in cooling water samples from conventional power plants as well as from nuclear power plants.

858 Professional Sample Processor – the universal sample changer

The 858 Professional Sample Processor is a flexible sample changer that is geared entirely to your requirements. It can be operated together with all intelligent Metrohm ion chromatographs.

The 858 Professional Sample Processor is equipped with a 786 Swing Head and a robotic transfer arm. This configuration makes it possible, firstly, to reach all positions on the sample rack. Secondly, it is also possible to reach external positions, which are used by the Liquid Handling Station. Furthermore, the 858 Professional Sample Processor can be combined with all Metrohm Inline Sample Preparation techniques. For example, your samples can be diluted and filtered directly during sample delivery.


Flexibility for a broad variety of uses

Thanks to the wide variety of equipment options and methods of working available, the 858 Professional Sample Processor is very flexible. Users can choose between a peristaltic pump or an 800 Dosino to introduce samples. Various injection methods, such as Full Loop, Internal Loop, Partial Loop or Pick-up, can be installed, and the sample can be transported into the injection valve in push or pull mode. The entire sample channel is metal-free and therefore also suitable for biological samples.

Sample volumes between 500 μ L and 500 mL can be processed. An optional injection valve directly on the sample changer provides for additional flexibility in sample prepa-

ration and introduction. Moreover, this optional injection valve enables using the 858 Professional Sample Processor as a fraction collector. In addition, a large selection of sample racks is available, so operators can use their own sample vessels. If a ready-made rack will not do, we will be happy to offer a customized rack. Up to 999 different sample positions can be accessed on Metrohm racks.

For Inline Sample Preparation or other tasks peripheral devices can be connected to and controlled by the sample changer. Optional peripheral devices include, for example, membrane pumps, magnetic stirrers, rod stirrers, remote connection boxes or additional 800 Dosinos.

The 858 Professional Sample Processor can be equipped with a variety of sample racks. Choose the rack that accommodates your amount of samples. You will find a selection of sample racks on page 20.

858 Professional Sample Processor

08

Intelligence matters

The 858 Professional Sample Processor is an integral part of Metrohm's intelligent ion chromatography. The intelligence of the various components minimizes operating errors and makes self-monitoring systems possible. Sample changer, sample rack and all connected peripheral devices are automatically recognized by the MagIC Net user software and are available directly with all functions without any further configuration.


More user comfort

The 858 Professional Sample Processor is characterized by exceptional ease of use. To increase its effectiveness, the 858 Professional Sample Processor can be used simultaneously with two or more analytical systems (multi-channel use). Thus you save laboratory space and reduce purchase costs.

Trace analysis as a special challenge

In trace and ultratrace analysis it is essential to eliminate any potential sources of contamination. With customized sample racks, users can place their sample vessels directly on the sample changer. Thus there is no need to transfer a sample from one vessel to another, which always involves a high risk of contamination. Impurities caused by an additional sample beaker are also ruled out.

The carryover from sample to sample can be minimized by various rinsing options. The rinsing of the sample needle from the inside and outside reduces carryover, for example, to $<0.001\%$ in the case of intelligent Partial Loop Injection. The Liquid Handling Station with integrated rinsing unit installed directly on the sample changer is used to this end. Alternatively, 300 mL beakers filled with ultrapure water can be defined as rinsing positions on the sample rack. Thus determinations in the mg/L range and trace determinations can be carried out one after the other in a sample series – without time-consuming rinsing procedures.


858 Professional Sample Processor in multi-channel use: An 858 Professional Sample Processor can introduce samples for two or more analytical systems.

Handwritten signature

Measuring pH value and conductivity in the sample stream

Measuring conductivity and the pH value in order to determine the dilution factor for a sample is a standard procedure in many laboratories. A single, fully automated IC system including both measurements can be achieved, if an 867 pH Module and an 856 Conductivity Module are combined with an 858 Professional Sample Processor. From the measured conductivity of the sample the soft-

ware calculates the required dilution factor. In combination with Metrohm Inline Dilution (MIDT), this procedure makes sure that the ionic concentration of the solution is ideal for the subsequent analysis by IC. If inline pH measurement is also included in this procedure, Metrohm enables users to determine the most important parameters in water analysis on a single system.

Integrated of inline determination of conductivity and pH combined with Metrohm Inline Dilution.

Handwritten signature


919 IC Autosampler plus – automation for medium-sized sample series

The 919 IC Autosampler plus is the ideal solution for simple automation tasks. Combined with a Metrohm intelligent ion chromatograph it enables sample transfer as well as frequently used Inline Sample Preparation steps such as Inline Ultrafiltration.

Generously equipped

56 positions for 11 mL respectively 2.5 mL sample vessels mean that the 919 IC Autosampler plus offers more than enough space. Due to its metal-free sample channel, both bio-chemical samples and samples with aggressive matrices can be reliably introduced. Sample transfer is provided either by a dual-channel peristaltic pump or by means of an 800 Dosino. Samples can be moved both in the push or pull mode. As a rinsing beaker can be easily installed on the rack, carryover effects from sample to sample are avoided.

Either simply used to introduce samples or with integrated Inline Ultrafiltration, the 919 IC Autosampler plus is an efficient solution for automation.

78x



Long-lived ceramic sample needle

The 919 IC Autosampler plus features a newly developed sample needle made from zirconium oxide and includes a tip made of PEEK material. The sample needle is rugged, hard to break and very durable. Moreover, minimizes carryover effects from sample to sample. The new sample needle is well-suited both for universal use in trace analysis as well as in highly contaminated samples.

Simple upgrade possible

If there is an increased number of samples to be analyzed, the 919 IC Autosampler can be easily upgraded at any point in time. In fact, it can be upgraded to feature the complete capacity and functionality of an 858 Professional Sample Processor: positions on the rack can be increased to 148, a Liquid Handling Station can be installed, and sample vessels of various shapes and sizes can be used. Moreover upgrading the 919 IC Autosampler plus enables the use of inline Dilution.

amw

The 919 IC Autosampler plus can be used with all Metrohm ion chromatographs.

Blas

863 Compact IC Autosampler – the smallest Metrohm sample changer


12

The 863 Compact IC Autosampler excels by its small footprint and impressive capability. The 863 Compact IC Autosampler automates the 940 Professional IC Vario, the 930 Compact IC Flex as well as the 883 Basic IC plus. The 863 Compact IC Autosampler, too, features a sample channel free of any metal. Thus the device is suitable both for trace analysis and for the analysis of aggressive main components. The 863 Compact IC Autosampler is the ideal tool for routine analysis.

Configuration for the small tasks in day-to-day laboratory work

The 863 Compact IC Autosampler is preconfigured and equipped with a sample rack for 36 positions. Its built-in peristaltic pump can be used to transfer the sample or for sample preparation steps. It can be operated at 7 speeds and in both directions; with this sample changer, too, sample can be supplied in the push and pull mode.

The sample changer can be controlled by all versions of MagIC Net or, alternatively, programmed easily by the integral keyboard. Moreover, you can automate several analytical systems operating independently of each other (multi-channel use). Even if you analyze only a limited number of samples, automation helps you save valuable time.



The 863 Compact IC Autosampler: small and reliable.

Very easy to use

The 863 Compact IC Autosampler is operated with the ion chromatography software MagIC Net. The software is available in a variety of languages; even users with basic skills will find it very easy to use.

«Remote» signals enable the 863 Compact IC Autosampler to be controlled even by Metrohm ion chromatographs that do not use MagIC Net or even by third-party instruments. The 863 Compact IC Autosampler is supplied with 4 predefined methods, which users can modify as they wish and store under a name of their choosing. To edit the methods, there is an integral keyboard combined with a large LCD screen. In the «remote» mode up to 99 samples can be processed in sequence. Here too, users can select their language, so operating errors are ruled out.

The 863 Compact IC Autosampler in combination with the 930 Compact IC Flex is the perfect system for routine analyses of small sample quantities.

889 IC Sample Center – the sample changer for small sample volumes

14

The 889 IC Sample Center is the perfect automation solution for small sample volumes. In combination with the optional cooling function it is the ideal system for sensitive biochemical samples.



The 889 IC Sample Center is a robust autosampler for high sample throughput and optimized for the challenges of the modern analytical laboratory. It works according to the x-y-z principle and with high-resolution injection control for precise sample delivery. A double-needle system enables vessel caps and septa to be pierced. A PEEK injection valve completes the automation system of the 889 IC Sample Center.

Fast and versatile

Less than 10 seconds from the start of the determination to injection – the 889 IC Sample Center is one of the fastest sample changers currently available in the market. Needless to say that a very high throughput is guaranteed.

The 889 IC Sample Center is also characterized by its versatility. Various injection methods can be selected: full-loop, partial-loop and pick-up.

The enormous flexibility of the system is also evident in the variability of the sample volume. The standard accessories include two 48-position-sample-racks for 300 and 700 μL sample vessels. Injection volumes from 1 to 100 μL can be selected. Other SBS standard sample racks that are common in liquid chromatography or microtiter plates can also be used with the 889 IC Sample Center. The system supports sample racks with 12 to 384 sample positions. Sample volumes may range from 1 μL to 10 mL.

The 889 IC Sample Center is the ideal sample changer for small sample volumes. The sample delivery is especially fast, which guarantees a high sample throughput.



758

Cool the sample for better precision

The 889 IC Sample Center is available with a cooling function. In this version a Peltier element ensures that samples are cooled to a stable temperature of minimum 4 °C. Due to the sophisticated cooling system no temperature gradient can develop; a constant temperature is

guaranteed. As a result, samples in which stability is critical are able to withstand long sample series unaffected; the precision of analytical results is increased. Thus the 889 IC Sample Center with cooling function is an attractive option not only for the biochemical and clinical sector.

Alufo

Just as all other Metrohm sample changers the 889 IC Sample Center can be combined with all intelligent Metrohm ion chromatographs. This offers high flexibility.

Alufo

814 USB Sample Processor – a versatile partner

16

The 814 USB Sample Processor extends the automation possibilities in ion chromatography. It guarantees a robust and highly precise introduction of samples and can be upgraded with various accessories. An 800 Dosino is used for sample delivery into the injection valve. The sample path of the 814 USB Sample Processor is also completely free of any metal and thus inert. Sensitive samples are protected.


For universal use

The 814 USB Sample Processor can be used not only for ion chromatography, but also for voltammetric and titrimetric applications. Various software packages such as *viva*, *tiamo*[™] and MagIC Net can be used to control the sample changer. Stand-alone operation with a touch panel is also possible. The great benefit for the user is that only one sample changer is needed for the different analytical methods.

Flexible equipment

The 814 USB Sample Processor is available with one work station and with two work stations. This opens up new possibilities for Inline Sample Preparation and for sample handling. For example an 814 USB Sample Processor with two working stations can introduce samples independently and without any contamination risk to a dual-channel system. The «Discover» function can be used to briefly remove the covers of large sample vessels (50, 75, and 250 mL) for sampling and seal them again.

Other advantages of the 814 USB Sample Processor are an optical beaker sensor, which virtually eliminates operating errors, an optional number of membrane pumps that can be used for rinsing procedures, and free choice of sample rack, allowing users to place their laboratory's specific sample vessels directly on the rack.

The 814 USB Sample Processor works through your samples continuously and with perfect reliability. An 800 Dosino is used to transfer samples, guaranteeing precise and contamination-free working.

815 Robotic USB Sample Processor XL – the largest Metrohm sample changer

The 815 Robotic USB Sample Processor XL is Metrohm's largest and most versatile sample changer for ion chromatography. Up to 228 samples, each of 11 mL, can be placed on the sample rack and processed in sequence. Your samples can be analyzed continuously not only overnight, but also, if necessary, over an entire weekend.


One sample changer for every application

Just like the 814 USB Sample Processor, the 815 Robotic USB Sample Processor XL can be used for voltammetry, titration and ion chromatography. It can be controlled by *viva*, *tiamo*™ and MagIC Net. All Metrohm sample racks can be used with the 815 Robotic USB Sample Processor XL. If a ready-made rack will not do, we will be happy to offer a customized solution. Hence, samples don't need to be transferred manually between different vials any longer.

Other advantages of the 815 Robotic USB Sample Processor XL are various rinsing options, to minimize the carryover from sample to sample, the «Discover» function, which allows even large sample volumes to be protected against dirt and contamination, and the beaker sensor, which virtually eliminates operating errors. The wide range of equipment available gives you the possibility to use one or two work stations and determine the number of membrane pumps yourself. Thus, the 815 Robotic USB Sample Processor XL offers enormous versatility for Inline Sample Preparation and Liquid Handling.

Cooling or heating of samples

The 815 Robotic USB Sample Processor XL can both cool (down to -10 °C) and heat (up to +55 °C) samples. By cooling, users can stabilize their samples and thereby achieve greater precision. In contrast, heating speeds up reactions; an effect that is used, for example, in sample preparation, to make unwanted components harmless.


The 815 Robotic USB Sample Processor XL with two work stations enables users to introduce samples and to collect different fractions of the analyzed samples. With 228 rack positions there is plenty of space available.

815 Robotic Soliprep for LC – automation for solid samples


18

The old cliché, according to which ion chromatography is only suitable for analyzing liquid samples, is no longer valid. Modern sample preparation and automation techniques also allow gases and solids to be analyzed. With the 815 Robotic Soliprep for LC, Metrohm offers a modern automation solution for solid samples.


The 815 Robotic Soliprep for LC is a version of the 815 Robotic USB Sample Processor XL, which is presented on page 17. With this system all manual steps can be automated in full.

Automatic sample preparation of various steps. In the first step, solvent is added, the sample comminuted and then homogenized. After that, sample is aspirated with a metal-free needle, filtered and finally transferred directly to an LC system by way of a connection port.

In the Soliprep system the sample rack can be adapted individually to the particular application. Different inserts provide the highest degree of flexibility.

The 815 Robotic Soliprep for LC – more than a sample changer


The 815 Robotic Soliprep for LC enables users to homogenize, extract, dilute and filter solid or liquid samples. Subsequently, the sample is transferred directly to an IC or HPLC system by way of a connector. Depending on the requirements of the application, the 815 Robotic Soliprep offers enormous flexibility in the configuration of the system and integrates all sample preparation steps into a single, continuous and fully automated analytical sequence. To this end there is a selection of sample racks available, which can be equipped with various inserts.

The equipment of the 815 Robotic Soliprep for LC
The 815 Robotic Soliprep for LC comprises not only an 815 Robotic USB Sample Processor XL with two work stations, but also the equipment for homogenizing the sample («Polytron» homogenizer and washing station for

the «Polytron»), for Liquid Handling (800 Dosinos and 772 Pump Unit) and for transferring the prepared sample to an LC system (connector with capillary). 24 solid samples in 120 mL vessels can be placed on the standard sample rack of the 815 Robotic Soliprep for LC. In addition, 30 sample vials (11 mL), 24 metal-free sample needles and 24 syringe filters are available there.

Applications from a variety of industrial sectors

With the 815 Robotic Soliprep for LC the scope of use of ion chromatography can also be extended to solid samples. Typical samples come from the pharmaceutical industry. However, this system is also ideally suited to the analysis of foods, feedstuffs and environmentally relevant samples such as sediments.

The 815 Robotic Soliprep for LC for Inline Sample Preparation of solid samples. The worked-up sample is transferred directly to an IC or HPLC system. This fully automatic analytical system minimizes manual steps and guarantees complete traceability.

Sample racks and sample vessels

A variety of sample racks can be used for Metrohm sample changers. A selection can be found in the table below. Plastic or glass sample vessels can be used together

with these racks. If your sample vessels are not listed, we would be happy to draw up a customized solution for you.

	Metrohm sample vessels															External sample vessels						
																Nalgene®			Eppendorf			
Material	Polypropylene										PE		Glass			PP, PE, FEP, PFA			PP			
Number of Sample positions on the rack	0.3 mL vials 6.2743.110	0.7 mL vials 6.2743.100	2.5 mL vials 6.2743.040	11 mL vials 6.2743.050	50 mL vessels 6.2747.010	75 mL vessels 6.1459.400	120 mL vessels 6.1459.300	200 mL vessels 6.1453.220	200 mL vessels 6.1459.310	250 mL vessels 6.1453.250	50 mL vessels 6.1608.100	300 mL vessels 6.1608.080	6 mL vessels 6.2419.000	75 mL vessels 6.1432.210	100 mL vessels 6.1608.050	250 mL vessels 6.1432.320	60 mL vessels	125 mL vessels	250 mL vessels	1 mL vials	1.5 mL vials	2 mL vials
Rack																						
6.2041.200	48	48																				
6.2041.210													12									
6.2041.310								12		12							12			12		
6.2041.340						24								24								
6.2041.350						48								48								
6.2041.360									12													
6.2041.370									14													
6.2041.380														14					14			
6.2041.390																		16				
6.2041.410			14	14							14											
6.2041.430			27	27							27											
6.2041.440			48	48							48											
6.2041.450			56	56	156																	
6.2041.480													19							159	159	159
6.2041.510			56	56							56											
6.2041.750			36	36																		
6.2041.760			54	54							54											
6.2041.800					100									100								
6.2041.820								28		28							28					
6.2041.830									28													
6.2041.850							59				59											
6.2041.860			228	228								12										
6.2041.900*						54									54							
6.9920.115*			95	95																		
6.9920.191											35						35		35			
6.9920.242	82	82																				

For use with 858 Professional Sample Processor, 814 USB Sample Processor and 815 Robotic USB Sample Processor XL

For use with 889 IC Sample Center

For use with 815 Robotic USB Sample Processor XL


PP ... Polypropylene

PE ... Polyethylene


FEP ... Perfluoroethylene propylene copolymer

PFA ... Perfluoroalkoxy copolymer

* Sample rack only for use in combination with water bath 6.2840.000 and drip pan 6.2711.070

Technical information

21

	858	919	863	889	814	815
Sample rack and sample vials	Variable: ≤ 999 positions standard rack: 148 × 2.5 or 11 mL vials	Variable: ≤ 720 positions standard rack: 56 × 2.5 or 11 mL vials	Fixed rack with 36 positions for 2.5 or 11 mL vials	Variable: ≤ 384 positions standard rack: 2 × 48 × 0.3 or 0.7 mL vials	Variable: ≤ 999 positions standard rack: 148 × 2.5 or 11 mL vials	Variable: ≤ 999 positions standard rack: 228 × 2.5 or 11 mL vials
Sample volume	0.5–500 mL	0.5–11 mL	0.5–11 mL	0.001–10 mL	0.5–500 mL	0.5–500 mL
Sample delivery	Dual-channel peristaltic pump or 800 Dosino	Dual-channel peristaltic pump or 800 Dosino	Single-channel peristaltic pump	Microliter syringe	800 Dosino	800 Dosino
External position accessible	Yes	No	No	No	Yes	Yes
Liquid Handling Station	Yes	No	No	Yes	Yes	Yes
Injection valve, e.g. for sample preparation	Yes	No	No	Yes	No	No
Metrohm Inline Sample Preparation «MISP»	All techniques	Some techniques	Limited selection	Limited selection	Some techniques	Some techniques
Injection method	Full Loop, Internal Loop, Partial Loop, Pick-up	Full Loop, Internal Loop, Partial Loop, Pick-up	Full Loop, Internal Loop, Partial Loop	Full Loop, Internal Loop, Partial Loop, Pick-up	Full Loop, Internal Loop, Partial Loop, Pick-up	Full Loop, Internal Loop, Partial Loop, Pick-up
Can be combined with	940, 930, 883	940, 930, 883	940, 930, 883	940, 930, 883	940, 930, 883	940, 930, 883
Software control by means of	MagIC Net, viva	MagIC Net, viva	MagIC Net	MagIC Net	MagIC Net, tiamo™, viva	MagIC Net, tiamo™, viva
Multi-channel use	Yes	Yes	Yes	Yes	Yes	Yes
Connection to PC	USB	USB	USB	USB	USB	USB
Control of peripheral devices via MSB	3 ×	3 ×	No	No	3 ×	3 ×
ZrO ₂ sample needle	Yes	Yes	Yes	No	Yes	Yes

Handwritten signature

Handwritten signature

Ordering information

22

	858 Professional Sample Processor
2.858.0010	858 Professional Sample Processor
2.858.0020	858 Professional Sample Processor with dual-channel peristaltic pump
2.858.0030	858 Professional Sample Processor with dual-channel peristaltic pump and 6-port injection valve
	919 IC Autosampler plus
2.919.0020	919 Autosampler plus incl. sample rack and vessels
6.5867.000	Upgrade Kit for 919 IC Autosampler plus
	863 Compact IC Autosampler
2.863.0010	863 Compact IC Autosampler including sample rack and sample vials
	889 IC Sample Center
2.889.0010	889 IC Sample Center including sample rack and sample vials
2.889.0020	889 IC Sample Center with cooling function including sample rack and sample vials
	814 USB Sample Processor
2.814.0010	814 USB Sample Processor with one work station and one membrane pump
2.814.0020	814 USB Sample Processor with one work station and two membrane pumps
2.814.0030	814 USB Sample Processor with one work station and no membrane pump
2.814.0110	814 USB Sample Processor with two work stations and two membrane pumps
2.814.0120	814 USB Sample Processor with two work stations and four membrane pumps
2.814.0130	814 USB Sample Processor with two work stations and no membrane pump
	815 Robotic USB Sample Processor XL
2.815.0010	815 Robotic USB Sample Processor XL with one work station and one membrane pump
2.815.0020	815 Robotic USB Sample Processor XL with one work station and two membrane pumps
2.815.0030	815 Robotic USB Sample Processor XL with one work station and no membrane pump
2.815.0110	815 Robotic USB Sample Processor XL with two work stations and two membrane pumps
2.815.0120	815 Robotic USB Sample Processor XL with two work stations and four membrane pumps
2.815.0130	815 Robotic USB Sample Processor XL with two work stations and no membrane pump
2.815.4110	815 Robotic Soliprep for LC
	Accessories for 814 USB Sample Processor and 815 Robotic USB Sample Processor XL
2.786.0010	786 Swing Head with robotic transfer arm (left)
2.786.0020	786 Swing Head with robotic transfer arm (right)
2.786.0040	786 Swing Head
	Liquid Handling
2.772.0110	772 Pump Unit
2.800.0010	800 Dosino
2.843.0020	843 Pump Unit with two membrane pumps
2.843.0120	843 Pump Unit with two peristaltic pumps
6.3032.120	Dosing Unit 2 mL
6.3032.150	Dosing Unit 5 mL
6.3032.210	Dosing Unit 10 mL
6.3032.220	Dosing Unit 20 mL
6.3032.250	Dosing Unit 50 mL

23

- 6.5330.100 IC equipment: Inline Dialysis
- 6.5330.110 IC equipment: Inline Ultrafiltration
- 6.5330.120 IC equipment: Inline Dilution
- 6.5330.140 IC equipment: MiPCT
- 6.5330.160 IC equipment: MiPCT-ME
- 6.5330.170 IC equipment: MiPuT
- 6.5330.180 IC equipment: MiPT

- 6.2841.120 Liquid Handling Station, left
- 6.2841.130 Liquid Handling Station, right
- 6.5330.130 IC equipment: Liquid Handling Station, left

Stirrers

- 2.741.0010 741 Magnetic stirrer for Liquid Handling Station
- 2.801.0010 801 Magnetic stirrer for external position
- 2.802.0020 802 Propeller stirrer



Sample vessels and closures

- 6.2743.050 11 mL sample vial, 2'000 pieces (PP)
- 6.2743.057 11 mL sample vial, 200 pieces (PP)
- 6.2743.040 2.5 mL sample vial, 2'000 pieces (PP)
- 6.2743.047 2.5 mL sample vial, 200 pieces (PP)
- 6.2743.070 Perforated stoppers for 11 and 2.5 mL sample vial, 2'000 pieces
- 6.2743.077 Perforated stoppers for 11 and 2.5 mL sample vial, 200 pieces

- 6.2743.100 0.7 mL sample vial, 1'000 pieces, PP
- 6.2743.107 0.7 mL sample vial, 100 pieces, PP
- 6.2743.110 0.3 mL sample vial, 1'000 pieces, PP
- 6.2743.117 0.3 mL sample vial, 100 pieces, PP
- 6.2743.120 Septum for 0.7 and 0.3 mL sample vial, 1'000 pieces
- 6.2743.127 Septum for 0.7 and 0.3 mL sample vial, 100 pieces

- 6.2820.000 Aluminum foil for closing large-volume sample vessels

Order numbers and information on sample racks and other sample vessels can be found in the table on page 20.

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- work 125

START WITH DETERMINATION:

Time	Sample	Valve	Fill
0.00	861 Advanced Compact IC	Valve	Fill
3.00	861 Advanced Compact IC	Valve	Inject
	838 Advanced Sample Processor	- Move	sample
	838 Advanced Sample Processor	- Lift	work
	838 Advanced Sample Processor	- Pump	180

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

CALIBRATION OF COMPONENT AMMONIUM

Method: C4 250.mvw
Equation: $y = 1.38189 \cdot x$
RSD: 0.000000
Correlation coefficient: 0.999900


CALIBRATION OF COMPONENT ETA

Method: C4 250.mvw
Equation: $y = 35.651 \cdot x$
RSD: 0.000000
Correlation coefficient: 0.999900

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




Sample data

Ident 5 ppb NH3+5 ppb MEA
 Sample type Проба
 Determination start 2018-11-18 16:22:34 UTC+2
 Method NH3 + MEA 15.11.2018
 Operator

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
Sample data

Ident 4 ppb NH3+4 ppb MEA
Sample type Проба
Determination start 2018-11-18 15:33:26 UTC+2
Method NH3 + MEA 15.11.2018
Operator

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

Sample data

Ident sid_4ppbNa+4ppbNH3
 Sample type Проба
 Determination start 2018-12-01 11:45:09 UTC+2
 Method Na + NH3 без подкисляване 29.11.2018
 Operator

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